



ClearSteth®

The ClearSteth® Hosted Software Solution is designed with your needs in mind. With an intuitive interface, teleauscultation is as simple as connecting your chest piece to the PC, opening the application and connecting to the clinician to remotely share high quality heart, lung, and body sounds. ClearSteth® is available as a standalone application or as an integrated application in our eNcounter® telemedicine software applications. Versatile and secure, users on various networks can access encrypted auscultation connections with the simple press of a button.



AUSCULTATIONS

CLEARSTETH®

Recommended bandwidth is 128 Kbps (standalone) 1Mbps download/upload (video conference)

SSL encryption to the eNcounter Server AES 256-bit for audio over UDP or TCP

Bell, Diaphragm, and Extended filters

Windows® 10, 11 Pro

Intel® Core® i5 (speeds of +2.4 GHz)

Minimum 4 GB RAM (recommended)

10.5 MB of disk space

USB 2.0, 3.0

+1 (480) 922-0044

globalmed.com

15023 N 73rd Street
Scottsdale, AZ 85260

