

### What's New for Users

- V-222536: Minimum password length upon reset will be 15 characters (increased from 8)
- V-222432: Enforce the limit of 3 consecutive invalid logon attempts by a user during a 15-minute time period

### What's New in ClearSteth®

- ClearSteth, UDP port fields will not accept non-numeric characters, or values outside of 0-65535, and will advise the users of invalid entries, and won't save the port value in the ClearSteth configuration settings
- ClearSteth, within eNcounter 3 now has the ability for the user to scroll based on scale of window being utilized

### What's New for Administrators

- All CAT 1 remediation updates
  - V-222432 (details above)
  - V-222536 (details above)
  - Upgraded server components from .Net6, which has gone end of life to .Net 8
- License Manager: Disabled "Quantity" field, so there can be no duplicate licenses created
- License Manager: Will not allow licenses to be created with an expiration date in the past
- Server Components with Updated Version Numbers:
  - GlobalMed\_PasswordReset\_1.0.4
  - GlobalMed\_ChangePassword\_1.0.1
  - GlobalMed\_eNcounterServer\_2.0.9.2
  - GlobalMed\_eNcounterAdmin\_1.1.11.0
  - GlobalMed\_LicenseService\_1.0.6.1
  - GlobalMed\_LicenseManager\_1.0.10.1
  - GlobalMed\_eNcounterView\_2.0.3
  - GlobalMed\_TurnCredentialsScrubber\_1.0.2
  - GlobalMed\_TurnCredentials\_1.0.2
  - GlobalMed\_TURNServer\_WebRTC\_UDP\_1.0.7.0
  - GlobalMed\_TURNServer\_UDP\_1.0.7.0
  - GlobalMed\_TURNServer\_TCP\_1.0.7.0
  - GlobalMed\_VersionedTextService\_1.0.2

© 2002-2024 GlobalMed Holdings, LLC. D.B.A. GlobalMed® All rights reserved. All trademarks are the property of their respective holder.

Reproduction or distribution of this document electronically or in print, is strictly prohibited without prior written permission from GlobalMed.