

eNcounter® 3.6.9.5 and ClearSteth 1.3.15.0

Release Notes



What's New for eNcounter®

- DICOM Manager Backup Management: Enhanced DICOM Manager with configurable backup retention settings.
- Updates to ClearSteth.
 - Updated ClearSteth background process (ClearSteth v1.3.15.0)
- General bug fixes - ECG incompatibility with parenthesis in patient ID, and TotalVitals hardware update.
- eView® bug fixes.

What's New for Administrators

We've introduced more transparency and granular control for system administrators within the Admin Tool:

- Configurable DICOM socket timeout settings per institution.
- Configurable DICOM backup retention per institution.

