

eNcounter® 3.6.9.6, ClearSteth® v1.3.16.0, eView v2.0.9, AdminTool v1.1.16, eCloud v3.0.1

Release Notes

What's New for Users: eNcounter, ClearSteth, eView

- Improved focus order to provide a logical and predictable tab sequence.
- Enhanced visible focus indicators to clearly show current keyboard focus.
- Improved usability via standard keyboard commands.
- Color Contrast & Visual Accessibility
- Updated color combinations to meet minimum contrast requirements for text, icons, and interactive controls.
- Improved readability for users with low vision and color-vision deficiencies.
- Screen Reader & Assistive Technology Support
- Improved screen reader compatibility.
- Enhanced announcements for form fields, buttons, dialogs, and status messages.
- Improved overall consistency of screen reader behavior across supported workflows.
- General Accessibility Enhancements
- Improved consistency of accessibility behavior across the application.
- Continued alignment with WCAG 2.x Level A and AA success criteria and Revised Section 508 standards.
- TotalVitals® USB Bug Fix
- When switching between cables TotalVitals will now function correctly and transfer vitals data as expected.

What's New for Administrators

- Allowing numerical values to be added to institution names