



GlobalMed

TotalECG™



Designed with simplicity, portability and accuracy in mind, this electrocardiogram system empowers providers to extend their cardiac care beyond traditional methodologies. The TotalECG is a resting ECG device incorporating a wireless 12-Lead ECG system, Windows®-based CardioVu software using the Glasgow 12-Lead ECG analysis program, and a compact ECG acquisition device.



CARDIAC EXAMINATIONS

TOTAL ECG®

Glasgow 12-lead Analysis at 500 sps

No Bluetooth® drivers required

One-click connection process

Optional 128-bit encryption

Runs on Windows®7 and 10

Works up to 33 ft. away from PC

Multiple devices will not
interfere up to 9

No power cord or outlet needed

+1 (480) 922-0044

globalmed.com

15023 N 73rd Street
Scottsdale, AZ 85260

