



**Transportable Exam Backpack**

The Transportable Exam Backpack is a lightweight solution that can be carried hands-free allowing clinicians to provide telemedicine services for patients in remote locations. This solution is equipped with protective foam inserts with compartments for individual medical devices. The backpack's soft case and multiple storage compartments allow users to fit additional devices and disposable components, safely storing up to a 17" laptop.

#### MEDICAL USE CASES

---



sonography



cardiology



dermatology



auscultations



audiology



ophthalmology



pulmonology



intraoral

## TRANSPORTABLE EXAM BACKPACK

---

#### WEIGHT

13.5lbs (6.12 kg)

#### POWER

Integrated 4 Port Hub

#### DIMENSIONS

21.6" x 13.4" x 7.0"

#### MOBILITY

30L Lightweight  
Water-Resistant Backpack

---

+1 (480) 922-0044

[globalmed.com](http://globalmed.com)

15023 N 73rd Street  
Scottsdale, AZ 85260

